




zot Antibody

Product Code	CSB-PA334425LA01VEX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P38442
Immunogen	Recombinant Vibrio cholerae serotype O1 Zona occludens toxin protein (1-399AA)
Raised In	Rabbit
Species Reactivity	Vibrio cholerae serotype O1
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zona occludens toxin (Zonular occludens toxin) (Zot), zot
Immunogen Species	Vibrio cholerae serotype O1
Target Names	zot

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: zot antibody at 3.5µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 52 kDa

Observed band size: 52 kDa