



HA-33 Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA338424LC01CLQ |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P46084 |
| Immunogen | Recombinant Clostridium botulinum Main hemagglutinin component protein (2-286AA) |
| Raised In | Rabbit |
| Species Reactivity | Clostridium botulinum |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Main hemagglutinin component (HA 33 kDa subunit) (HA1), HA-33, antP-33 ha1 |
| Immunogen Species | Clostridium botulinum |
| Target Names | HA-33 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |