



HA-33 Antibody, Biotin conjugated

Product Code	CSB-PA338424LD01CLQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P46084
Immunogen	Recombinant Clostridium botulinum Main hemagglutinin component protein (2-286AA)
Raised In	Rabbit
Species Reactivity	Clostridium botulinum
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Main hemagglutinin component (HA 33 kDa subunit) (HA1), HA-33, antP-33 ha1
Immunogen Species	Clostridium botulinum
Target Names	HA-33
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.