



# botF Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA341192HC01CLQ
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P30996
<b>Immunogen</b>	Recombinant Clostridium botulinum Botulinum neurotoxin type F protein (1-436AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Clostridium botulinum
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Botulinum neurotoxin type F (BoNT/F) (EC 3.4.24.69) (Bontoxilysin-F) [Cleaved into: Botulinum neurotoxin F light chain; Botulinum neurotoxin F heavy chain], botF
<b>Immunogen Species</b>	Clostridium botulinum
<b>Target Names</b>	botF
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.