



VP3 Antibody, Biotin conjugated

Product Code	CSB-PA346512LD01CIF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P54095
Immunogen	Recombinant Chicken anemia virus Apoptin protein (1-120AA)
Raised In	Rabbit
Species Reactivity	Chicken anemia virus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Apoptin, VP3
Immunogen Species	Chicken anemia virus
Target Names	VP3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.