



ompA Antibody, Biotin conjugated

Product Code	CSB-PA356881LD01DSG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06597
Immunogen	Recombinant Chlamydia trachomatis serovar L2 Major outer membrane porin protein (23-394AA)
Raised In	Rabbit
Species Reactivity	Chlamydia trachomatis serovar L2
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Major outer membrane porin (MOMP), ompA, omp1L2
Immunogen Species	Chlamydia trachomatis serovar L2
Research Area	Microbiology
Target Names	ompA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.