




vpX Antibody

Product Code	CSB-PA356976LA01HLZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06939
Immunogen	Recombinant Human immunodeficiency virus type 2 subtype A Protein Vpx protein (1-112AA)
Raised In	Rabbit
Species Reactivity	Human immunodeficiency virus type 2 subtype A
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein Vpx (Viral protein X) (X ORF protein), vpx
Immunogen Species	Human immunodeficiency virus type 2 subtype A
Target Names	vpx

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: vpx at 2.6µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 32 kDa

Observed band size: 32 kDa