




NS1 Antibody

Product Code	CSB-PA357083LA01HPM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07298
Immunogen	Recombinant Human parvovirus B19 Non-capsid protein NS-1 protein (1-255AA)
Raised In	Rabbit
Species Reactivity	Human parvovirus B19
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Non-capsid protein NS-1 (NCVP1) (Non-structural protein NS1), NS1
Immunogen Species	Human parvovirus B19
Target Names	NS1

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: NS1 antibody at 2.7µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 45 kDa

Observed band size: 45 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.