



NS1 Antibody, HRP conjugated

Product Code	CSB-PA357083LB01HPM
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07298
Immunogen	Recombinant Human parvovirus B19 Non-capsid protein NS-1 protein (1-255AA)
Raised In	Rabbit
Species Reactivity	Human parvovirus B19
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Non-capsid protein NS-1 (NCVP1) (Non-structural protein NS1), NS1
Immunogen Species	Human parvovirus B19
Target Names	NS1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.