




pmpD Antibody

Product Code	CSB-PA530746LA01DSB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O84818
Immunogen	Recombinant Chlamydia trachomatis Probable outer membrane protein PmpD protein (1244-1531AA)
Raised In	Rabbit
Species Reactivity	Chlamydia trachomatis
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable outer membrane protein PmpD (Polymorphic membrane protein D), pmpD
Immunogen Species	Chlamydia trachomatis
Target Names	pmpD

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: pmpD antibody at 3.9µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 47 kDa

Observed band size: 47 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.