



pmpD Antibody, HRP conjugated

Product Code	CSB-PA530746LB01DSB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O84818
Immunogen	Recombinant Chlamydia trachomatis Probable outer membrane protein PmpD protein (1244-1531AA)
Raised In	Rabbit
Species Reactivity	Chlamydia trachomatis
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable outer membrane protein PmpD (Polymorphic membrane protein D), pmpD
Immunogen Species	Chlamydia trachomatis
Target Names	pmpD