



pmpD Antibody, Biotin conjugated

Product Code	CSB-PA530746LD01DSB
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O84818
Immunogen	Recombinant Chlamydia trachomatis Probable outer membrane protein PmpD protein (1244-1531AA)
Raised In	Rabbit
Species Reactivity	Chlamydia trachomatis
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable outer membrane protein PmpD (Polymorphic membrane protein D), pmpD
Immunogen Species	Chlamydia trachomatis
Target Names	pmpD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.