




# NS Antibody

<b>Product Code</b>	CSB-PA604972LA01IER
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q0HD54
<b>Immunogen</b>	Recombinant Influenza A virus Non-structural protein 1 protein (1-237AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Influenza A virus
<b>Tested Applications</b>	ELISA, WB; Recommended dilution: WB:1:500-1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Non-structural protein 1 (NS1) (NS1A), NS
<b>Immunogen Species</b>	Influenza A virus
<b>Target Names</b>	NS

## Image



### Western Blot

Positive WB detected in Recombinant protein

All lanes: NS antibody at 5µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 32 kDa

Observed band size: 32 kDa

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.