




P Antibody

Product Code	CSB-PA768100LA01VBH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86132
Immunogen	Recombinant Vesicular stomatitis Indiana virus Protein C' protein (1-67AA)
Raised In	Rabbit
Species Reactivity	Vesicular stomatitis Indiana virus
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein C', P
Immunogen Species	Vesicular stomatitis Indiana virus
Target Names	P

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: P antibody at 2.7µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 14 kDa

Observed band size: 14 kDa