



P Antibody, Biotin conjugated

Product Code	CSB-PA768100LD01VBH
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86132
Immunogen	Recombinant Vesicular stomatitis Indiana virus Protein C' protein (1-67AA)
Raised In	Rabbit
Species Reactivity	Vesicular stomatitis Indiana virus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein C', P
Immunogen Species	Vesicular stomatitis Indiana virus
Target Names	P
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.