




vcp Antibody

Product Code	CSB-PA772050LA01DIL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7ZU99
Immunogen	Recombinant Zebrafish Transitional endoplasmic reticulum ATPase protein (704-806AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transitional endoplasmic reticulum ATPase (TER ATPase) (EC 3.6.4.6) (Protein CDC48) (Valosin-containing protein) (VCP), vcp, cdc48
Immunogen Species	Zebrafish
Research Area	Neuroscience; Metabolism; Signal transduction
Target Names	vcp

Image



Western Blot

Positive WB detected in: Zebrafish tissue, Zebrafish brain tissue, Zebrafish heart tissue, Zebrafish eye tissue

All lanes: vcp antibody at 1.8µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 90 kDa

Observed band size: 90 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.