



# L Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA806333LD01VBE
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q8B0H0
<b>Immunogen</b>	Recombinant Vesicular stomatitis Indiana virus RNA-directed RNA polymerase L protein (598-784AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Vesicular stomatitis Indiana virus
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	RNA-directed RNA polymerase L (Protein L) (Large structural protein) (Replicase) (Transcriptase) [Includes: RNA-directed RNA polymerase (EC 2.7.7.48); mRNA (guanine-N(7)-)-methyltransferase (EC 2.1.1.56); GDP polyribonucleotidyltransferase (EC 2.7.7.88); Cap-specific mRNA (nucleoside-2'-O-)-methyltransferase 2 (EC 2.1.1.296)], L
<b>Immunogen Species</b>	Vesicular stomatitis Indiana virus
<b>Target Names</b>	L