



L Antibody, Biotin conjugated

Product Code	CSB-PA806333LD01VBE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8B0H0
Immunogen	Recombinant Vesicular stomatitis Indiana virus RNA-directed RNA polymerase L protein (598-784AA)
Raised In	Rabbit
Species Reactivity	Vesicular stomatitis Indiana virus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RNA-directed RNA polymerase L (Protein L) (Large structural protein) (Replicase) (Transcriptase) [Includes: RNA-directed RNA polymerase (EC 2.7.7.48); mRNA (guanine-N(7)-)-methyltransferase (EC 2.1.1.56); GDP polyribonucleotidyltransferase (EC 2.7.7.88); Cap-specific mRNA (nucleoside-2'-O-)-methyltransferase 2 (EC 2.1.1.296)], L
Immunogen Species	Vesicular stomatitis Indiana virus
Target Names	L
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.