



VP1 Antibody, FITC conjugated

Product Code	CSB-PA859915LC01CID
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99153
Immunogen	Recombinant Chicken anemia virus Capsid protein (288-369AA)
Raised In	Rabbit
Species Reactivity	Chicken anemia virus
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Capsid protein (CA1) (Coat protein), VP1
Immunogen Species	Chicken anemia virus
Target Names	VP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.