




PG_1424 Antibody

Product Code	CSB-PA886515LA01EYA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9RQJ2
Immunogen	Recombinant Porphyromonas gingivalis Peptidylarginine deiminase protein (44-556AA)
Raised In	Rabbit
Species Reactivity	Porphyromonas gingivalis
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peptidylarginine deiminase (EC 3.5.3.-)
Immunogen Species	Porphyromonas gingivalis
Target Names	PG_1424

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: PG_1424 antibody at 3.7µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 63 kDa

Observed band size: 63 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.