




Recombinant Variola virus A27L protein (A27L)

Product Code	CSB-EP2678VAR
Abbreviation	Recombinant Variola virus A27L protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q89694
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Variola virus (isolate Human/India/Ind3/1967) (VARV) (Smallpox virus)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MEKEIQRLRDRIMDLDRQLNECKRNGNGTSSEEVNRLKTRISDLKQSLEICKDK KSELYSAYKTELGRVREQISNLQESLRERESDKTDSYYRRELTRERNKIVELE KELNKCDFVNHVKYIDEINSKKNRISDLERQLTACKSNGCGNGDMDQYKREIE SLKRELTKCRRGSNGSHSDCEYYDEEARDCVKS
Research Area	Others
Source	E.coli
Target Names	A27L
Expression Region	1-194aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	27.0 kDa
Protein Length	Full Length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.