



TARDBP Antibody, Biotin conjugated

Product Code	CSB-PA023129ED01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13148
Immunogen	Recombinant Human TAR DNA-binding protein 43 protein (1-396AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	TAR DNA-binding protein 43 (TDP-43), TARDBP, TDP43
Immunogen Species	Homo sapiens (Human)
Research Area	Microbiology
Target Names	TARDBP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.