



NAIP Antibody, Biotin conjugated

Product Code	CSB-PA615662LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13075
Immunogen	Recombinant Human Baculoviral IAP repeat-containing protein 1 protein (147-242AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Baculoviral IAP repeat-containing protein 1 (Neuronal apoptosis inhibitory protein), NAIP, BIRC1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	NAIP