



NAIP Antibody, Biotin conjugated

Product Code	CSB-PA615662LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13075
Immunogen	Recombinant Human Baculoviral IAP repeat-containing protein 1 protein (147-242AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Baculoviral IAP repeat-containing protein 1 (Neuronal apoptosis inhibitory protein), NAIP, BIRC1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	NAIP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.