




Recombinant Human A disintegrin and metalloproteinase with thrombospondin motifs 2 (ADAMTS2), partial

Product Code	CSB-EP001308HU
Abbreviation	Recombinant Human ADAMTS2 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O95450
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	RRRARRHAADDDYNIEVLLGVDDSVVQFHHGKEHVQKYLTLNMNIVNEIYHDES LGAHINVVLVRIILLSYGKSMLEIGNPSQSLENVCRWAYLQQKPDTGHDEYH DHAIFLTRQDFGPSGMQGYAPVTGMCHPVRCTLNHEDGFSSAFVVAHETGH VLGMEHDGQGNRCGDEVRLGSIMAPLVQAAFHRFHWSRCSQQELSRYLHSY DCLLDDPFAHDWPALPQLPGLHYSMNEQC
Research Area	Cancer
Source	E.coli
Target Names	ADAMTS2
Protein Names	Recommended name: A disintegrin and metalloproteinase with thrombospondin motifs 2 Short name= ADAM-TS 2 Short name= ADAM-TS2 Short name= ADAMTS-2 EC= 3.4.24.14 Alternative name(s): Procollagen I N-proteinase Short
Expression Region	254-492aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	32.2 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.