




Recombinant Human Developmentally-regulated GTP-binding protein 2 (DRG2)

Product Code	CSB-EP2817HU
Abbreviation	Recombinant Human DRG2 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A8MZF9
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MGILEKISEIEKEIARTQKNKATEYHLGLLKAKLAKYRAQLLEPSKSASSKGEFG DVMKSGDARVALIGFPSVGKSTFLSLMTSTASEAASYEFTTLTCIPGVIEYKGA NIQLLDLPGIIEGAAQGGKGRGRQVIAVARTADVIIMMLDATKGEVQRSLLLEKELE SVGIRLNKHKPNIFYKPKKGGGISFNSTVTLTQCSEKLVQLILHEYKIFNAEVLFR EDCSPDEFIDVIVGNRVYMPCLYVYNKIDQISMEEVDRLARKPNSVVISCGMKL NLDYLLEMLWEYLALTCIYTKKRGQRPDFTDAILRKGASVEHVCHRIHRSLAS QFKAPAPSTVRSGWA
Research Area	Cell Biology
Source	E.coli
Target Names	DRG2
Expression Region	1-343aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	44.1 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.