




# Recombinant Human Proto-oncogene tyrosine-protein kinase Src (SRC)

<b>Product Code</b>	CSB-EP022650HU
<b>Abbreviation</b>	Recombinant Human SRC protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P12931
<b>Storage Buffer</b>	Tris-based buffer,50% glycerol
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥ 85% as determined by SDS-PAGE.
<b>Sequence</b>	GSNKS PKDASQRRRSLEPAENVHGAGGGAFPASQTPSKPASADGHRGPSA AFAPAAAEPKLFGGFNSSDVTVTSPQRAGPLAGGVTTFVALYDYESRTETDLSF KKGERLQIVNNTGDWWLAHSLSTGQTGYIPSNYVAPSDSIQAEWYFGKITR RESERLLLNAENPRGTFVLVRESETTKGAYCLSVSDFDNAKGLNVKHYKIRKLD SGGFYITSRTQFNSLQQLVAYYSKHADGLCHRLTTCPTSKPQTQGLAKDAW EIPRESLRLEVKLGQGCFCGEVWMGTWNGTTRVAIKTLKPGTMSPEAFLQEAQ VMKCLRHEKLVQLYAVVSEEPYIVTEYMSKGSLLDFLKGETGKYLRLPQLVDM AAQIASGMAYVERMNYVHRDLRAANILVGENLVCKVADFGFLARLIEDNEY TAR QGAKFPIKWTAPEAALYGRFTIKSDVWSFGILLTELTTKGRVPYPGMVNREVL D QVERGYRMPCPPECPESLHDLMCQCWRKEPEERPTFEYLQAFLEDYFTSTEP QYQPGENL
<b>Research Area</b>	Cancer
<b>Source</b>	E.coli
<b>Target Names</b>	SRC
<b>Protein Names</b>	Recommended name: Proto-oncogene tyrosine-protein kinase Src EC= 2.7.10.2 Alternative name(s): Proto-oncogene c-Src pp60c-src Short name= p60-Src
<b>Expression Region</b>	2-536aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 10xHis-tagged
<b>Mol. Weight</b>	64.2 kDa
<b>Protein Length</b>	Full Length of Mature Protein
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.