




# Recombinant Streptococcus pyogenes serotype M28 Holo-[acyl-carrier-protein] synthase (acpS)

<b>Product Code</b>	CSB-EP669769SBAF
<b>Abbreviation</b>	Recombinant Streptococcus pyogenes serotype M28 acpS protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q48RM7
<b>Storage Buffer</b>	Tris-based buffer,50% glycerol
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Streptococcus pyogenes serotype M28 (strain MGAS6180)
<b>Purity</b>	Greater than 85% as determined by SDS-PAGE.
<b>Sequence</b>	MIVGHGIDLQEISAIEKVYQRNPRFAQKILTEQELAI FESFPYKRRLSYLAGRWA GKEAFAKAIGTGIGRLTFQDIEILNDVRGCPILTKSPFKGNSFISISHSGNYVQAS VILEDKK
<b>Research Area</b>	Others
<b>Source</b>	E.coli
<b>Target Names</b>	acpS
<b>Protein Names</b>	Recommended name: Holo-[acyl-carrier-protein] synthase Short name= Holo-ACP synthase EC= 2.7.8.7Alternative name(s): 4'-phosphopantetheinyl transferase AcpS
<b>Expression Region</b>	1-118aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 10xHis-tagged and C-terminal Myc-tagged
<b>Mol. Weight</b>	18.3 kDa
<b>Protein Length</b>	Full Length
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.