





Recombinant Human Meteorin (METRN)

Product Code	CSB-EP892322HU
Abbreviation	Recombinant Human METRN protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UJH8
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	<p>GYSEERC¹SWRGSGLTQEPGSVGLALACAEGAVEWLYPAGALRLTLGGPDP RARPGIACLRPVRPFAGAQVFAERAGGALELLLAEGPGPAGGRCV²⁰RWGP RRALFLQATPHQDISRRVAAFRFELREDGRPELPPQAHGLGVDGACR³⁰PCSDA ELLAACTSDFVIHGIIHGVTHDVELQESVITVVAARVLRQTPLPFQAGR⁴⁰SGDQ GLTSIRTP⁵⁰LRCGVHPGPGTFLFMGWSRFGEARLGCAPRFQEFRRAYEA⁶⁰ARAA HLHPCEVALH</p>
Research Area	Neuroscience
Source	E.coli
Target Names	METRN
Protein Names	Recommended name: Meteorin
Expression Region	24-293aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	36.5 kDa
Protein Length	Full Length of Mature Protein


Image



Based on the SEQUEST from database of E.coli host and target protein, the LC-MS/MS Analysis result of CSB-EP892322HU could indicate that this peptide derived from E.coli-expressed Homo sapiens (Human) METRN.



Based on the SEQUEST from database of E.coli host and target protein, the LC-MS/MS Analysis result of CSB-EP892322HU could indicate that this peptide derived from E.coli-expressed Homo sapiens (Human) METRN.



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.