




Recombinant Human Meteorin (METRN)

Product Code	CSB-EP892322HUa0
Abbreviation	Recombinant Human METRN protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9UJH8
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	GYSEERC ^{SW} RGSGLTQEPG ^{SV} Q ^L ALACAEGAVEWLYPAGALRLTLGGPDP RARPGIACLRPVRPFAGAQVFAERAGGALELLLAEGPGPAGGRCVRWGP RRALFLQATPHQDISRRVAAFRFELREDGRPELPPQAHGLGVDGACRPCSDA ELLLAACTSD ^{FVI} HGIIHG ^{VTH} DVELQESVITVVAARVLRQT ^{PPL} FQAGRSGDQ GLTSIRTP ^{LR} CGVHPGPGTFLFMGWSRFGEARLGCAPRFQEFRRAYEAARAA HLHPCEVALH
Research Area	Neuroscience
Source	E.coli
Target Names	METRN
Expression Region	24-293aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	33.1 kDa
Protein Length	Full Length of Mature Protein

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.