




Recombinant Human adenovirus C serotype 5 DNA-binding protein (DBP), partial

Product Code	CSB-EP365892HIL
Abbreviation	Recombinant Human adenovirus C serotype 5 DBP protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P03265
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Human adenovirus C serotype 5 (HAdV-5) (Human adenovirus 5)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	SVPIVSAWEKGMEAAARALMDKYHVDNDLKANFKLLPDQVEALAAVCKTWLNE EHRGLQLTFTSKKTFVTMMGRFLQAYLQSFAEVTYKHHEPTGCALWLHRC AEI EGELKCLHG SIMINKEHVIEMDVTSENGQRALKEQSSKAKIVKNRWGRNVVQI SNTDARCCVHDAAC PANQFSGKSCGMFFSEGAKAQVAFKQIKAFMQALYPNA QTGHGHLLMPLRCECNSKPGHAPFLGRQLPKLTPFALSNAEDLDADLISDKSV LASVHHPALIVFQCCNPVYRNSRAQGGGPNCDFKISAPDLLNALVMVRSLWSE NFTELPRMVVPEFKWSTKHQYRNVSLPVAHSDARQNPFD F
Research Area	Cancer
Source	E.coli
Target Names	DBP
Protein Names	Recommended name: Early E2A DNA-binding protein
Expression Region	174-529aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	43.9 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.