



Recombinant Rat Serotransferrin (Tf)

Product Code	CSB-YP023412RA
Abbreviation	Recombinant Rat Serotransferrin protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P12346
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	VPDKTVKWC AVSEHENTK CISFRDHMKT VLPADG PRLACVKKTSYQDCIKAIS GGEADAITLDGGWVYDAGLTPNNLKPVA AEFYGSLEHPQTHYLAVAVVKKGT DFQLNQLQGKKSCHTGLGRSAGWIIPIGILLFCNLPEPRKPLEKAVASFFSGSCV PCADPVAFPQLCQLCPGCGCSPTQPFPGYVGA FKCLRDGGGDVAFVKHTTIF EVL PQKADR DQYELLCLDNTRK PVDQYEDCYLARIPSHAVVARNGDGKEDLIW EILKVAQEHFGKGKSKDFQLFGSPLGKDLLFKDSAFGLLRVPPRMDYRLYLGH SYVTAIRNQREGVCPEGSIDSAPVKWCALSHQERAKCDEWSVSSNGQIECES AESTEDCIDKIVNGEADAMSLDGGHAYIAGQCGLVPVMAENYDISSCTNPQSD VFPKGYAVAVVKASDSSINWNNLKGKKSCHTGVDRTAGWNIPMGLLFSRINH CKFDEFFSQGCAPGYKKNSTLCDLCIGPAKCAPNNREGYNGYTGA FQCLVEK GDVAFVKHQTVLENTNGKNTAAWAKDLKQEDFQLLCPDGT KKPVTEFATCHL AQAPNHVVVSRKEKAARVSTVLTAQKDLFWKGD KCTGNFCLFRSSTKDLLF RDDTKCLTKLPEGTTYEEYLGAEY LQAVGNIRK CSTSRLLEACTFHKS
Research Area	Epigenetics and Nuclear Signaling
Source	Yeast
Target Names	Tf
Protein Names	Recommended name: Serotransferrin Short name= Transferrin Alternative name(s): Beta-1 metal-binding globulin Liver regeneration-related protein LRRG03 Siderophilin
Expression Region	20-698aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	76
Protein Length	Full Length of Mature Protein
Image	

