




Recombinant Vesicular stomatitis Indiana virus Matrix protein (M)

Product Code	CSB-YP356390VBG
Abbreviation	Recombinant Vesicular stomatitis Indiana virus Matrix protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P04876
Storage Buffer	Tris-based buffer, 50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Vesicular stomatitis Indiana virus (strain Glasgow) (VSIV)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MSSLKKILGLKGKGKKSKKLGIAPPPYEEDTSMEYAPSAPIDKSYFGVDEMDTH DPNQLRYEKSFFTVKMTVRSNRPFRITYSDVAAAVSHWDHMYIGMAGKRPFYK ILAFLGSSNLKATPAVLADQGQPEYHAHCEGRAYLPHRMGKTPPMLNVPEHF RRPFNIGLYKGTIELTMTIYDDESLEAAPMIWDHFNSSKFSDFREKALMFGLIVE EEASGAWVLDSVRHSKWASLASSF
Research Area	Signal Transduction
Source	Yeast
Target Names	M
Protein Names	Recommended name: Matrix protein
Expression Region	1-237aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	28.8
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.