




Recombinant Human Dihydropyrimidinase-related protein 3 (DPYSL3)

Product Code	CSB-EP622756HUa0
Abbreviation	Recombinant Human DPYSL3 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q14195
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MSYQGKKNIPRITSDRLLIKGGRIVNDDQSFYADIYMEDGLIKQIGDNLIVPGGV KTIEANGKMVIPGGIDVHTHFQMPYKGMTTVDDFFQGTKAALAGGTTMIIDHVV PEPESSLTEAYEKWREWADGKSCCDYALHVDITHWNDSVKQEVQNLIKDKGV NSFMVYMAYKDLYQVSNTELYEIFTCLGELGAIAQVHAENGDIIAQEQTMRMLEM GITGPEGHVLSRPEELEAEAVFRAITIASQTNCPLYVTKVMSKSAADLISQARKK GNVVFGEPIASLGIDGTHYWSKNWAKAAAFVTSPPLSPDPTTPDYINSLLASG DLQLSGSAHCTFSTAQKAIGKDNFTAPEGTNGVEERMSVIWDKAVATGKMDE NQFVAVTSTNAAKIFNLYPRKGRISVGSDSLVIWDPDAVKIVSAKNHQSAAEY NIFEGMELRGAPLVVICQGKIMLEDGNLHVTQGAGRIFPCSPFSDYVYKRIKAR RKMADLHAVPRGMYPDGPVFDLTTTPKGGTPAGSARGSPTRPNPPVRNLHQS GFSLSGTQVDEGVRSASKRIVAPPGGRSNITSLS
Research Area	Neuroscience
Source	E.coli
Target Names	DPYSL3
Expression Region	1-570aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	66.0 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.