





Recombinant Mouse Ventricular zone-expressed PH domain-containing protein 1 (Veph1), partial

Product Code	CSB-EP025844MO
Abbreviation	Recombinant Mouse Veph1 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A1A535
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	AMTESTFPQQKDLEQLQLHLEEVRRFFDVFGFSETAGAWQCFMCNNPEKATVV NQDGQPLIEGKLKEKQVRWKFIKRWKTHYFTLAGNQLLFQKGKSKDDPDSP ELSKVQSVKAVAKKRRDRSLPRAFEIFTDSKTYVFKAKDEKNAAEELWQCIN VALAQAKERESREVTYYL
Research Area	Others
Source	E.coli
Target Names	Veph1
Protein Names	Recommended name: Ventricular zone-expressed PH domain-containing protein 1 Alternative name(s): Protein melted homolog
Expression Region	660-833aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	27.7 kDa
Protein Length	Partial


Image



Based on the SEQUEST from database of E.coli host and target protein, the LC-MS/MS Analysis result of CSB-EP025844MO could indicate that this peptide derived from E.coli-expressed Mus musculus (Mouse) Veph1.



Based on the SEQUEST from database of E.coli host and target protein, the LC-MS/MS Analysis result of CSB-EP025844MO could indicate that this peptide derived from E.coli-expressed Mus musculus (Mouse) Veph1.



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.