




Recombinant Rat Probable E3 ubiquitin-protein ligase makorin-3 (Mkrn3)

Product Code	CSB-BP736743RA
Abbreviation	Recombinant Rat Mkrn3 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	XP_001057262.2
Product Type	Recombinant Proteins
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MEESTAPMEAYEAAGAEAGASSSSSSPPQFEASGASAVASSAPLQEASGLA PFLMAPVPANRRAASLRPAPAEGGGARSGPERNTGSWTKQILCRYYLHGQC KEGDNCRYSHDLSGRRKARGGQDSQPRASADRGPKMATHWEPPTQEVAEA PPTASSSSLPLIGSAAERGFSEAEIDNAGIGSAAERGFPEAEIDNAGLAAGAAG GAGAEGWEGAIEFVPGQPYPYRGRMIPPHGPEAPLQSPEIEREHMAMGMGMPL PLCRYAARGQCLRGDRCAYPHGEICDMCGQQALHPWDAAQQEAHRRACVE AHERDMELSFVQRSMKVCVICMEVVYEKADPSDRRFILFSCNHTYCLKCI RRWRSATQFENRISKSCPQCRVSSGFVIPSEFWVEEEEEKEKLVQQYKEGMS QKACRYFAGGLGHCPFGFCFYKHEYPEGWRDQPPRPDGGGSSSAYWHQV LEPVQLREGSVLFKSRKKEHSVLRRLANQLLKLKLLCLRGSFSISDDRWLLLQYQL EEYFNLNL
Research Area	Neuroscience
Source	Baculovirus
Target Names	Mkrn3
Expression Region	60.0 kDa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Mol. Weight	60.0 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.