





Recombinant Human adenovirus C serotype 5 DNA-binding protein (DBP), partial

Product Code	CSB-BP365892HIL
Abbreviation	Recombinant Human adenovirus C serotype 5 DBP protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P03265
Product Type	Recombinant Proteins
Immunogen Species	Human adenovirus C serotype 5 (HAdV-5) (Human adenovirus 5)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	SVPIVSAWEKGMEAAARALMDKYHVDNDLKANFKLLPDQVEALAAVCKTWLNE EHRGLQLTFTSKKTFVTMMGRFLQAYLQSF AEVTYKHHEPTGCALWLHRC AEI EGELKCLHGSIMINKEHVIEMDVTSENGQRALKEQSSKAKIVKNRWGRNVVQI SNTDARCCVHDAACPANQFSGKSCGMFFSEGAKAQVAFKQIKAFMQALYPNA QTGHGHLLMPLRCECNSKPGHAPFLGRQLPKLTPFALSNAEDLDADLISDKSV LASVHHPALIVFQCCNPVYRNSRAQGGGPNCDFKISAPDLLNALVMVRSWSE NFTELPRMVVPEFKWSTKHQYRNVSLPVAHSDARQNPFD
Research Area	Cancer
Source	Baculovirus
Target Names	DBP
Protein Names	Recommended name: Early E2A DNA-binding protein
Expression Region	174-529aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Mol. Weight	42.3 kDa
Protein Length	Partial


Image



Based on the SEQUEST from database of Baculovirus host and target protein, the LC-MS/MS Analysis result of CSB-BP365892HIL could indicate that this peptide derived from Baculovirus-expressed Human adenovirus C serotype 5 (HAdV-5) (Human adenovirus 5) DBP.



Based on the SEQUEST from database of Baculovirus host and target protein, the LC-MS/MS Analysis result of CSB-BP365892HIL could indicate that this peptide derived from Baculovirus-expressed Human adenovirus C serotype 5 (HAdV-5) (Human adenovirus 5) DBP.



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.