




# Recombinant Mouse Lipoprotein lipase (Lpl)

<b>Product Code</b>	CSB-BP013065MO
<b>Abbreviation</b>	Recombinant Mouse Lpl protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P11152
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	≥ 85% as determined by SDS-PAGE.
<b>Sequence</b>	ADAGRDFSDIESKFALRTPEDTAEDTCHLIPGLADSVSNCHFHNHSSKTFVVIHG WTVTGMYESWVPKLVAALYKREPDSNVIVVDWLYRAQQHYPVSAGYTKLVG NDVARFINWMEEEFNYPDLNVHLLGYSLGAAHAAGVAGSLTNKKVNTRITGLDPA GPNFEYAEAPSRLSPDDADFVDVLHTFTRGSPGRSIGIQKPVGHVDIYPNGGT FQPGCNIGEAIRVIAERGLGDVDQLVKCSHERSIHLFIDSLLEENPSKAYRCN SKEAFEKGLCLSCRKNRCNNLGYEINKVRAKRSSKMYLKTRSQMPYKVFHYQ VKIHFSGTEDGKQHNQAFEISLYGTVAESENIPFTLPEVSTNKTYSFLLIYTEVDIG ELLMMLKLVWSDSYFSWPDWWSSPSFVIERIRVKAGETQKKVIFCAREKVS HLQKGD SAVFVKCHDKSLKKSG
<b>Research Area</b>	Cancer
<b>Source</b>	Baculovirus
<b>Target Names</b>	Lpl
<b>Protein Names</b>	Recommended name: Lipoprotein lipase Short name= LPL EC= 3.1.1.34
<b>Expression Region</b>	52.8 kDa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 10xHis-tagged
<b>Mol. Weight</b>	28-474aa
<b>Protein Length</b>	Full Length of Mature Protein
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.