




# Recombinant Human Diacylglycerol kinase alpha (DGKA)


<b>Product Code</b>	CSB-BP006832HU
<b>Abbreviation</b>	Recombinant Human DGKA protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P23743
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥ 85% as determined by SDS-PAGE.
<b>Sequence</b>	MAKERGLISPSDFAQLQKYMEYSTKKVSDVLKLFEDGEMAKYVQGDAIGYEGF QQFLKIYLEVDNVPRHLSLALFQSFETGHCLNETNVTKDVVCLNDVSCYFSLLE GGRPEDKLEFTFKLYDTRNGILDSSEVDKILQMMRVAEYLDWDVSELRPILQ EMMKEIDYDGGSGSVSQAEWVRAGATTVPLLVLGLEM TLKDDGQHMWRPKR FPRPVYCNLCCESSIGLGKQGLSCNLCKYTVHDQCAMKALPCEVSTYAKSRKDI GVQSHVWVRGGCESGRCDRCQKKIRIYHSLTGLHCVWCHLEIHDDCLQAVG HECDCGLLRDHILPPSSIYPSVLASGPDRKNSKTSQKTMDLNLSTSEALRIDP VPNTHPLLVFVNPKSGGKQGRVLWKFQYILNPRQVFNLLKDGPEIGLRLFKD VPDSRILVCGGDGTWILETIDKANLPVLPVAVLPLGTGNDLARCLRWGGG YEGQNLAKILKDLEMSKVHMDRWSVEVIPQQTEEKSDPVPFQIINNYFSIGVD ASIAHRFHIMREKYPEKFNSRMKNKLWYFEFATSE SIFSTCKKLEESLTVEICGK PLDLSNLSLEGIAVLNIPSMHGGSNLWGDTRRPHGDIYGINQALGATAKVITDP DILKTCVPDLSDKRLEVVGLEGA IEMGQIYTKLKNAGRRLAKCSEITFHHTKTLP MQIDGEPWMQTPCTIKITHKNQMPMLMGPPPRSTNFFGFLS
<b>Research Area</b>	Signal Transduction
<b>Source</b>	Baculovirus
<b>Target Names</b>	DGKA
<b>Protein Names</b>	Recommended name: Diacylglycerol kinase alpha Short name= DAG kinase alpha EC= 2.7.1.107 Alternative name(s): 80 kDa diacylglycerol kinase Diglyceride kinase alpha Short name= DGK-alpha
<b>Expression Region</b>	86.7 kDa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 10xHis-tagged and C-terminal Myc-tagged
<b>Mol. Weight</b>	86.7 kDa
<b>Protein Length</b>	Full Length




## Image



Based on the SEQUEST from database of Baculovirus host and target protein, the LC-MS/MS Analysis result of CSB-BP006832HU could indicate that this peptide derived from Baculovirus-expressed Homo sapiens (Human) DGKA.



Based on the SEQUEST from database of Baculovirus host and target protein, the LC-MS/MS Analysis result of CSB-BP006832HU could indicate that this peptide derived from Baculovirus-expressed Homo sapiens (Human) DGKA.



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

## Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.