




Recombinant Human E3 ubiquitin-protein ligase pellino homolog 1 (PELI1)

Product Code	CSB-BP856932HUb0
Abbreviation	Recombinant Human PELI1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96FA3
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MFSPDQENHPSKAPVKYGELIVLGYNGSLPNGDRGRRKSRFALFKRPKANGV KPSTVHIACTPQAAKAISNKDQHSISYTLRAQTVVVEYTHDSNTDMFQIGRST ESPIDFVVTDTVPGSQSNSDTQSVQSTISRACRIICERNPPFTARIYAAGFDSS KNIFLGEKAAKWKTSDGQMDGLTTNGVLVMHPRNGFTEDSKPGIWREISVCG NVFSLRETRSAQQRGKMVEIETNQLQDGSLIDLGCATLLWRTAEGLSHTPTVK HLEALRQEINAARPQCPVGFNTLAFPSMKRKDVVDEKQPWVYLNCGHVHGYH NWDGNKEERDGDRECPMCRSVGPYVPLWLWGCEAGFYVDAGPPTHAFSPCG HVCSEKTTAYWSQIPLPHGTHTFHAACPFCAHQLAGEQGYIRLIFQGPLD
Research Area	Signal Transduction
Source	Baculovirus
Target Names	PELI1
Expression Region	1-418aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Mol. Weight	48.8 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.