




Recombinant Human Complement C1q tumor necrosis factor-related protein 3 (C1QTNF3)

Product Code	CSB-BP883621HU
Abbreviation	Recombinant Human C1QTNF3 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9BXJ4
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	QDEYMESPQTGGLPPDCSKCCHGDYSFRGYQGPPGPPGPPGIPGNHGNNG NNGATGHEGAKGEKGDKGD LGPRGERGQHGPKEGKGYPGIPPELQIAFMAS LATHFSNQNSGIIFSSVETNIGNFFDVM TGRFGAPVSGVYFFTF SMMKHEDVE EVYVYLMHNGNTVFSMYSYEMKGKSDTSSNHAVLKLAKGDEVWLRMGNGAL HGDHQRFFSTFAGFLLFETK
Research Area	Metabolism
Source	Baculovirus
Target Names	C1QTNF3
Protein Names	Recommended name: Complement C1q tumor necrosis factor-related protein 3 Alternative name(s): Secretory protein CORS26
Expression Region	28.1 kDa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	28.1 kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.