





Recombinant Human Actin-like protein 8 (ACTL8)

Product Code	CSB-YP887982HU
Abbreviation	Recombinant Human ACTL8 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9H568
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MAARTVIIDHGSGFLKAGTAGWNEPQMVFNPVNYLPCKENPGPSYARRRVSL GIDICHDPDTFSYPIERGRILNWEQVQYLWSFVLENHRREQVPPVIITETPLREP ADRRKMLEILFELLHVPSVLLADQLQMSLYASGLLTGVVVDVSGYGLTRVQPFH QGRPLPASGKTLEFAGQDLSAYLLKSLFKEDCDRRCLEFQLETVAVTQMKNKCY VPQNLGEALDFRERQQSALDESNTYQLPDGSRVELTPMQRVAPEMFFSPQVF EQPGPSIPRAIVESVESCEISLRPLLVSHVMACGGNTLYPGFTKRLFRELMGDH VSSTKATVWEGSNRNFVWLGASVVAHLSTYQSEWMSREEYGEHMRM
Research Area	Signal Transduction
Source	Yeast
Target Names	ACTL8
Protein Names	Recommended name: Actin-like protein 8 Alternative name(s): Cancer/testis antigen 57 Short name= CT57
Expression Region	1-366aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Mol. Weight	43.9 kDa
Protein Length	Full Length

Image



Based on the SEQUEST from database of Yeast host and target protein, the LC-MS/MS Analysis result of CSB-YP887982HU could indicate that this peptide derived from Yeast-expressed Homo sapiens (Human) ACTL8.



Based on the SEQUEST from database of Yeast host and target protein, the LC-MS/MS Analysis result of CSB-YP887982HU could indicate that this peptide derived from Yeast-expressed Homo sapiens (Human) ACTL8.



Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.