




Recombinant Human Protein FAM217A (FAM217A)

Product Code	CSB-EP809009HU
Abbreviation	Recombinant Human FAM217A protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8IXS0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MGRRNENCANSLRVSNISQENLSHWNLDSEVPVSENKNLPAGRDGAAGGKIN KNYLEIPVEQLMLEPNLSVHSQKSTQNSKQGIFQLWNCPLNEGSTIEKREFKK SSVETGFNVINHPIRVFTLNHPLTIASVDKQVGPYPGLPMLGLCWPYADGDFF KNRNEIHVSSCSTIENNDGETLPAPNWNLKHGNSSVEENFTDESDLSENEKTN DTLLSYFKKVDLNLKPETIKNVEEPPFTEEPNEVFPYPDFLPPPFSAALDLHNLALS KSDNWKVTVDPAETSVEHLITRLELERLQHMTIQKERPRLQTTFFCTPAVTERP SSSKATPKVRQPKLCDLSLQIPCVDKSQEKSKNNSGCKLEQNALKRNWSN AGKYRWNSRPLSLKSSSTPKQLIETYDKNPKSSILSPCQELSFKPTIGHTNQSM VKMVSTRCLPWRSPMPVSPILTFPENQKEEIKAPKRNFGTKKKLYRQNIVLN RPFSIQKLNCLSPSLIAKDKCCSPIEQK
Research Area	Cell Biology
Source	E.coli
Target Names	FAM217A
Protein Names	Recommended name: Protein FAM217A
Expression Region	1-508aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	61.4 kDa
Protein Length	Full Length

Image



Based on the SEQUEST from database of E.coli host and target protein, the LC-MS/MS Analysis result of CSB-EP809009HU could indicate that this peptide derived from E.coli-expressed Homo sapiens (Human) FAM217A.

