





Recombinant Treponema phagedenis Treponemal membrane protein B (tmpB)

Product Code	CSB-EP333469THS
Abbreviation	Recombinant Treponema phagedenis tmpB protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P29720
Product Type	Recombinant Protein
Immunogen Species	Treponema phagedenis
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	ISYDDNEYSRKSRAYTQLAEKAYDAGEYDTAIEYSQLAENFAQESAHEYIKRMM ARNEAEDAMNKARTRYAWAKQQKADKNYPTEYLIAGEAIKAGGIAFDNKNYDV AVTCAEKALESLKTVEPEDKVIAKAAADKAAAEEKAAKEKAAREKSAKDKAKEK AAKEKAAKDKAKEKAAKEKAAKDKAKEKAAKEKAAREMAAKEKAAKDCAA KEEAARKAAEEAAARKAAEEAAARKAAEEEAARIAEEEAARKAAEEEAARKA AEEALYNEKGKVLVLPSEYKVLTWKLDRECFWNIKNPAVYNDPFWWRKLYEA NKDKIPESNNPDWVEAETILVIPSIRGERREGLYDPDVKYQALPKR
Research Area	Others
Source	E.coli
Target Names	tmpB
Protein Names	Recommended name: Treponemal membrane protein B Alternative name(s): Antigen tmpB
Expression Region	22-384aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	47.4 kDa
Protein Length	Full Length of Mature Protein

Image





Based on the SEQUEST from database of E.coli host and target protein, the LC-MS/MS Analysis result of CSB-EP333469THS could indicate that this peptide derived from E.coli-expressed *Treponema phagedenis* tmpB.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.