




Recombinant Mycobacterium bovis EspC protein homolog (BQ2027_MB3645C)

Product Code	CSB-EP353261MVH
Abbreviation	Recombinant Mycobacterium bovis BQ2027_MB3645C protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P65088
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Mycobacterium bovis (strain ATCC BAA-935 / AF2122/97)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MTENLTVQPERLGV LASHHDNAAVDASSGVEAAAGLGESVAITHGPYCSQFN DTLNVYLTAHNALGSSLHTAGVDLAKSLRIA AKIYSEADEAWRKAIDGLFT
Research Area	Others
Source	E.coli
Target Names	BQ2027_MB3645C
Protein Names	Recommended name: EspC protein homolog
Expression Region	1-103aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Mol. Weight	16.3 kDa
Protein Length	Full Length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.