



Recombinant Human RB1-inducible coiled-coil protein 1 (RB1CC1), partial

Product Code	CSB-EP855061HU1
Abbreviation	Recombinant Human RB1CC1 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8TDY2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	AIQTALKEFKLEREVVEKELLEKVKHLENQIAKSPAIDSTRGDSSSLVAELQEKL QEEKAKFLEEQLEEQEKRKNEEMQNVRTSLIAEQQTNFNTVLTREKMRKENIIN DLSDKLKSTMQQQERDKDLIESLSEDRARLLEEKKKLEEEVSKLRSSSFVPSP YVATAPELYGACAPELPGESDRSAVETADEGRVDSAMETSMMSVQENIHMLS EEKQRIMLLERTLQLKEEENKRLNQRLMSQSMSSVSSRHSEKIAIRDFQVGDL VLIILDERHDNYVLFTVSPTLYFLHSESLPALDLKPGEGASGASRRPWVLGKVM EKEYCQAKKAQNRFKVPLGTFYRVKAVSWNKKV
Research Area	Others
Source	E.coli
Target Names	RB1CC1
Expression Region	1241-1594aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	47.7 kDa
Protein Length	Partial


Image



Based on the SEQUEST from database of E.coli host and target protein, the LC-MS/MS Analysis result of CSB-EP85061HU1 could indicate that this peptide derived from E.coli-expressed Homo sapiens (Human) RB1CC1.



Based on the SEQUEST from database of E.coli host and target protein, the LC-MS/MS Analysis result of CSB-EP85061HU1 could indicate that this peptide derived from E.coli-expressed Homo sapiens (Human) RB1CC1.



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.