




# Recombinant Human Ecto-NOX disulfide-thiol exchanger 2 (ENOX2), partial

<b>Product Code</b>	CSB-EP007677HU2a0
<b>Abbreviation</b>	Recombinant Human ENOX2 protein, partial
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q16206
<b>Storage Buffer</b>	Tris-based buffer,50% glycerol
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥ 85% as determined by SDS-PAGE.
<b>Sequence</b>	GILIQFEQIVAVYHSASKQKAWDHFTKAQRKNISVWCKQAEIIRNIHNDELMI RREEEMEMSDDEIEEMTETKETEEESALVSQAEALKEENDSLRWQLDAYRNEV ELLKQEQGKVHREDDPNKEQQLKLLQQALQGMQQHLLKVQEEYKKKEAELEK LKDDKLQVEKMLENLKEKESCASRLCASNQDSEYPLEKTMNSSPIKSEREALL VGIISTFLHVHPFGASIEYICSYLHRLDNKICTSDVECLMGRQLQHTFKQEMTGVG ASLEKRWKFCGFEGFLKLT
<b>Research Area</b>	Cancer
<b>Source</b>	E.coli
<b>Target Names</b>	ENOX2
<b>Expression Region</b>	327-610aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-tagged
<b>Mol. Weight</b>	37.1 kDa
<b>Protein Length</b>	Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.