




Recombinant Mouse Dedicator of cytokinesis protein 8 (Dock8), partial

Product Code	CSB-EP807451MO
Abbreviation	Recombinant Mouse Dock8 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8C147
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus (Mouse)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	KSYQASPDRLRLTWLQNMAEKHTKKKCFTEAMCLVHAAALVAEYLSMLEDHS YLPVGSVSFQNISSNVLEESAVSDDTLSPDEDGVCSGRYFTESGLVGLLEQAA ELFSTGGLYETVNEVYKLVIPILEAHRDFRKLSTHDKLQKAFDNIINKDHKRMF GTYFRVGFYGSRFGDLDEQEFVYKEPAITKLPEISHRLEGFYGQCFGAEFVEVI KDSTPVDKTKLDPNKAYIQITFVEPYFDEYEMKDRVITYFEKNFNLRRFMYTTPF TLEGRPRGELHEQHRRTVLTTMHAFPIKTRIRVSQKEEFVLTPIEVAIEDMK KKTLQLAVATHQEPPDAKMLQMVLQGSVGATVNQGPLEVAQVFLAEIPADPKL YRHHNKLRLCFKEFIMRCGEAVEKNRRLITAEQREYQQELKKNYNKLRLDSLRLP MIERKIP
Research Area	Others
Source	E.coli
Target Names	Dock8
Protein Names	Recommended name: Dedicator of cytokinesis protein 8
Expression Region	1633-2067aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	55.8 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.