





Recombinant Mouse Methylcytosine dioxygenase TET2 (Tet2), partial

Product Code	CSB-EP678395MO
Abbreviation	Recombinant Mouse Tet2 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q4JK59
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	RISLVLYRHKNLFLPKHCLALWEAKMAEKARKEEECGKNGSDHVSQKNHGKQ EKREPTGPQEPSYLRFIQSLAENTGSVTTDSTVTTSPYAFTQVTGPYNTFV
Research Area	Epigenetics and Nuclear Signaling
Source	E.coli
Target Names	Tet2
Protein Names	Recommended name: Methylcytosine dioxygenase TET2 EC= 1.14.11.n2 Alternative name(s): Protein Ayu17-449
Expression Region	1810-1912aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-KSI-tagged
Mol. Weight	27.0 kDa
Protein Length	Partial


Image



Based on the SEQUEST from database of E.coli host and target protein, the LC-MS/MS Analysis result of CSB-EP678395MO could indicate that this peptide derived from E.coli-expressed Mus musculus (Mouse) Tet2.



Based on the SEQUEST from database of E.coli host and target protein, the LC-MS/MS Analysis result of CSB-EP678395MO could indicate that this peptide derived from E.coli-expressed Mus musculus (Mouse) Tet2.



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.