



Recombinant Mouse Nuclear receptor subfamily 2 group F member 6 (Nr2f6)

Product Code	CSB-EP016058MO
Abbreviation	Recombinant Mouse Nr2f6 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P43136
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus (Mouse)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MAMVTGGWGDPPGGDTNGVDKAGGSYPRATEDDSASPPGATSDAEPGDEER PGLQVDCVVC GDKSSGKH YGVFTCEGCKSFFKRSIRRNLSYTCRSNRDCQID QHHRNQ CQYCR LKKCFRVGM RKEAVQRGRIPHALPGPAACSPPGATGVEPF TGPPVSELIAQLLRAEPYPAAGRFGGGGAVLGIDNVCELAARLLFSTVEWARH APFFPELPAADQVALLRLSWSELFVLNAAQAALPLHTAPLLAAAGLHAAPMAA ERAVAFMDQVRAFQE QVDKLGRLQVDAAEYGCLKAIALFTP DACGLSDPAHV ESLQEKAQVALTEYVRAQYPSQPQRFRGRLRLPALRAVPASLISQLFFMRLV GKTPIETLIRDMLLSGSTFNWPYGSG
Research Area	Epigenetics and Nuclear Signaling
Source	E.coli
Target Names	Nr2f6
Protein Names	Recommended name: Nuclear receptor subfamily 2 group F member 6 Alternative name(s): COUP transcription factor 3 Short name= COUP-TF3 V-erbA-related protein 2 Short name= EAR-2
Expression Region	1-390aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	49.4 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.